

Why Is My DNA or RNA Extraction Yield Low?


Common Causes and Solutions

Introduction

High-quality DNA and RNA are essential for PCR, RT-qPCR, sequencing, cloning, and other molecular biology applications.

However, researchers frequently encounter low extraction yields, poor nucleic acid purity, degradation, and inconsistent results.

The underlying causes may involve sample handling, incomplete lysis, inefficient purification, or unsuitable extraction conditions.



1. Low DNA or RNA Yield

The Problem

The final nucleic acid concentration is lower than expected.

Possible Causes

- Insufficient starting material.
- Incomplete cell or tissue lysis.
- Inefficient nucleic acid binding to the purification matrix.
- Sample loss during washing or transfer.
- Incomplete elution.
- Extraction chemistry unsuitable for the sample type.

Practical Solutions

Verify the recommended sample input and extraction protocol.

Ensure that the lysis procedure is appropriate for the biological material and that purification reagents are prepared correctly.

Check the manufacturer's recommended elution conditions and volume.

Practical tip: A smaller elution volume may increase the measured nucleic acid concentration without increasing the total recovered amount.

2. DNA or RNA Degradation

The Problem

Extracted nucleic acids appear fragmented or show reduced integrity.

Possible Causes

- Nuclease contamination.
- Delayed sample processing.
- Inappropriate storage conditions.
- Repeated freeze-thaw cycles.
- Excessive mechanical disruption.

Practical Solutions

Use appropriate nuclease-free materials and validated sample-preservation procedures.

Minimize unnecessary handling and follow sample-specific storage recommendations.

For RNA, take additional precautions to minimize RNase exposure.

3. Poor Nucleic Acid Purity

The Problem

Purity measurements suggest contamination, or downstream reactions perform poorly.

Possible Causes

- Residual proteins.
- Phenol or other extraction reagents.
- Guanidinium salts.
- Ethanol carryover.
- Incomplete washing.
- Contaminating DNA in RNA preparations.

Practical Solutions

Review the purification and washing steps. Ensure that residual wash buffers are removed according to the kit instructions.

Use appropriate purity and integrity assessments before downstream applications.

Important: Absorbance ratios alone cannot confirm that a nucleic acid sample is free from all contaminants or suitable for every downstream assay.

4. Poor PCR or RT-qPCR Performance After Extraction

The Problem

Extracted nucleic acids are measurable, but amplification fails or produces inconsistent results.

Possible Causes

- Residual PCR inhibitors.
- Degraded template.
- Inaccurate concentration measurements.
- DNA contamination in RNA samples.
- Insufficient template quality for the assay.

Practical Solutions

Assess sample integrity and consider appropriate template dilution to investigate inhibition.

For RNA-based experiments, evaluate the need for DNase treatment and include suitable controls for genomic DNA contamination.

5. Inconsistent Extraction Results Between Samples

The Problem

Similar samples produce substantially different yields or purity measurements.

Possible Causes

- Variation in sample quantity or quality.
- Inconsistent lysis efficiency.
- Pipetting errors.
- Differences in processing time.
- Sample-specific inhibitor content.

Practical Solutions

Standardize sample collection, input quantity, storage, and processing conditions.

Include appropriate extraction controls and document deviations from the validated protocol.

Quick Troubleshooting Reference

Problem	Possible Cause	Recommended Check
Low yield	Incomplete lysis or inefficient recovery	Review sample input and lysis
Degraded nucleic acid	Nuclease activity or poor storage	Evaluate handling and integrity
Poor purity	Reagent carryover	Review washing and purification
PCR inhibition	Residual contaminants	Assess inhibitors and template dilution
Variable results	Inconsistent sample processing	Standardize extraction workflow

Conclusion

Successful DNA and RNA extraction depends on appropriate sample preparation, effective lysis, efficient purification, and careful nucleic acid handling.

A systematic evaluation of yield, purity, integrity, and downstream performance can help identify extraction problems and improve experimental reliability.

Relevant Laboratory Products

EMDA International supplies DNA and RNA extraction kits, spin columns, nuclease-free consumables, and molecular biology reagents.

Email Now: sales@emdainternational.com

For research use only. Follow manufacturer instructions and sample-specific laboratory protocols.